



CFT-308 Infrared Thermometer

The Choicemmed infrared Thermometer is a reusable thermometer intended for the intermittent determination of human body temperature on people of all ages.

It can be used by consumers in the household environment and by healthcare provider.

FEATURES

- Quick measurement, less than 1 second
- Easy operation with intelligence.
- Multi-functional, It not only supports for human, but is able to take room/object 20 groups of memories, easy to recall.
- Age temperature distinction.
- Safe and hygienic . No touch helps minimize spreading of germs, Completely safe for use on Children and adults.
- Fever indicating system, the thermometer for fever will indicate temperature status with corresponding face expression and beeps. Get better noticed of your temperature, urge you to take corresponding actions.
- Auto shut-down and power-saving.
- CFDA approval

TECHNICAL SPECIFICATIONS

Application Range

Adult and child (ear test mode only for the child above 3 months).

Measurement

Measurement range:	32.0° C ~ 42.9° C (89.6° F ~ 109.2° F)	Measurement method	No Contact-Type
--------------------	--	--------------------	-----------------

Precision:	±0.2°C/0.4°F	Resolution:	0.1°C/°F
------------	--------------	-------------	----------

Display

Type:	OLED	Dimension:	151X35X50.2MM
-------	------	------------	---------------

Parameters:

Forehead temperature	Weight:	About 103 g (without battery)
----------------------	---------	-------------------------------

Power Requirements

Battery:	AAAX2
----------	-------

Data Transmitting

Blue Tooth 4.0	Waterproof
----------------	------------

Environment Requirements

Operation Temperature:	16° C ~ 35° C	Operation Humidity:	15%~85% no condensation
------------------------	---------------	---------------------	-------------------------

Storage Temperature:	- 20° ~ 55° C	Storage Humidity:	15%~85%, no condensation
----------------------	---------------	-------------------	--------------------------

Accessories

AAAx2 Batteries	Color box	One user manual
-----------------	-----------	-----------------