

U80M

Upper Arm Electronic Blood Pressure Monitor With Bluetooth



Irregular heartbeat (IHB) indicator



Certified according to EU product regulations

Specification

Product Specification	
Measuring Accuracy	Pressure: ± 3 mmHg Pulse: $\pm 5\%$ of reading
Measuring Range	Pressure: 0 to 299 mmHg Pulse: 40 to 199 pulses/min
Inflation	Automatic by micro air pump
Deflation	Automatic electronic control valve
Automatic Power-Off	Approx. 3 minutes
Operating Environment	+5°C to +40°C 15% - 85%RH
Storage Environment	-20°C to +55°C 10% - 85%RH

Application / Features



Year/Month/Date /Time function



2x90 sets memory



Large LCD




3 times average




Automatic power off

Easy to use




Step 1-preparation

Wrap the cuff around the upper arm of your left hand.



Step 2-measuring

Press Start/Stop button to start the measuring.



Memory check

Press MEM button to view the previous measurement result.

CE ROHS ISO BSCI FDA